### Important Safety Information

- Read the full manual before use. Safety may be impaired if it is not followed.
- Only connect to the voltage and frequency shown on the serial number plate. An earth is required.
- Always ensure the bath is filled with a liquid rated for the working temperature. See manual for liquid selection.
- Ensure the T100/TC120 is positioned so that the plug is easily accessible for disconnection.
- Before moving or cleaning, allow the unit to cool down and disconnect from the power supply.
- Safety can only be assured if the T100/TC120 is securely installed on a suitable vessel. See manual for guidance.


### Features

1. Alarm Light
2. Display
3. F (function) Button
4. Pump Outlet Plate (*TC120 only*)
5. Heater On Light
6. Over-Temperature Dial (*TC120 only*)
7. S (select) Button
8. Dial
9. Float: low liquid level protection

### Setting the temperature

1. Press S: Display will flash
2. Rotate dial to desired temp
3. Press S to store value


### Setting the over-temp (TC120 only)

Rotate dial to desired temp cut-out (above set point)

### Using the external pump (TC120 only)

Ensure the TC120 is switched off. Connect hoses securely to the pump outlet plate as shown left.

**Always use pump connectors and hoses that are suitable for the operating temperature and liquid used.**

**Never disconnect any pipes or hoses while they contain very hot or very cold liquids or while the TC120 is pumping.**

### Safety cut out and alarm light

The T100/TC120 is fitted with a float switch which monitors the liquid level during operation. Should the liquid level drop below a certain level, the safety cut out will be activated, and alarm light will appear. Top up the bath with liquid until at least a minimum level to resume operation.

### Cleaning

Clean the unit with a damp cloth, using only water. Do not use chemical cleaning agents.

Consult manual for full instructions.

### Service & repair

For service or repair please contact our Service Department or your local distributor. All returned units must be accompanied by a Return Materials Authorisation (RMA) number, obtainable by contacting the Grant service department.

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