


Overall dimensions




Packed weight: 4kg
Product weight: 3kg

Packed dimensions;
A 445mm
B 325mm
C 300mm

Product dimensions;
A 335mm
B 215mm
C 200mm

Working Area



D	298mm (at opening)
D	281mm (at tray height)
E	148mm (at opening)
E	131mm (at tray height)
F	150mm
G	132mm
H	150mm
I	191mm

Vessel capacity with tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	8
Duran 21-801: 150ml	62mm	110mm	8
Duran 21-801: 250ml	70mm	143mm	6s
Duran 21-801: 500ml	86mm	181mm	0
Duran 21-801: 750ml	95mm	203mm	0
Duran 21-801: 1000ml	101mm	230mm	0
Duran 21-801: 2000ml	136mm	265mm	0

Vessel capacity without tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	8
Duran 21-801: 150ml	62mm	110mm	8
Duran 21-801: 250ml	70mm	143mm	6s
Duran 21-801: 500ml	86mm	181mm	3
Duran 21-801: 750ml	95mm	203mm	2
Duran 21-801: 1000ml	101mm	230mm	2
Duran 21-801: 2000ml	136mm	265mm	0

s - staggered to fit maximum number