


Overall dimensions




Packed weight: 7.5kg
Product weight: 6kg

Packed dimensions;
A 470mm
B 475mm
C 325mm

Product dimensions;
A 360mm
B 380mm
C 225mm

Working Area



D	323mm (at opening)
D	306mm (at tray height)
E	298mm (at opening)
E	281mm (at tray height)
F	150mm
G	132mm
H	150mm
I	241mm

Vessel capacity with tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	23s
Duran 21-801: 150ml	62mm	110mm	16
Duran 21-801: 250ml	70mm	143mm	14s
Duran 21-801: 500ml	86mm	181mm	8s
Duran 21-801: 750ml	95mm	203mm	3s
Duran 21-801: 1000ml	101mm	230mm	2
Duran 21-801: 2000ml	136mm	265mm	0

Vessel capacity without tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	23s
Duran 21-801: 150ml	62mm	110mm	16
Duran 21-801: 250ml	70mm	143mm	14s
Duran 21-801: 500ml	86mm	181mm	9
Duran 21-801: 750ml	95mm	203mm	8s
Duran 21-801: 1000ml	101mm	230mm	5s
Duran 21-801: 2000ml	136mm	265mm	4

s - staggered to fit maximum number