


Overall dimensions




Packed weight: 11kg
Product weight: 9kg

Packed dimensions;
A 445mm
B 700mm
C 375mm

Product dimensions;
A 335mm
B 590mm
C 275mm

Working Area



D	298mm (at opening)
D	281mm (at tray height)
E	501mm (at opening)
E	485mm (at tray height)
F	150mm
G	132mm
H	150mm
I	234mm

Vessel capacity with tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	37s
Duran 21-801: 150ml	62mm	110mm	28
Duran 21-801: 250ml	70mm	143mm	24s
Duran 21-801: 500ml	86mm	181mm	10
Duran 21-801: 750ml	95mm	203mm	5s
Duran 21-801: 1000ml	101mm	230mm	0
Duran 21-801: 2000ml	136mm	265mm	0

Vessel capacity without tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	37s
Duran 21-801: 150ml	62mm	110mm	28
Duran 21-801: 250ml	70mm	143mm	24s
Duran 21-801: 500ml	86mm	181mm	15
Duran 21-801: 750ml	95mm	203mm	12s
Duran 21-801: 1000ml	101mm	230mm	12s
Duran 21-801: 2000ml	136mm	265mm	6

s - staggered to fit maximum number