


Overall dimensions




Packed weight: 16.5kg
Product weight: 14.5kg

Packed dimensions;
A 455mm
B 880mm
C 470mm

Product dimensions;
A 335mm
B 770mm
C 370mm

Working Area



D	302mm (at opening)
D	281mm (at tray height)
E	635mm (at opening)
E	635mm (at tray height)
F	191mm
G	180mm
H	191mm
I	270mm

Vessel capacity with tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	50s
Duran 21-801: 150ml	62mm	110mm	45s
Duran 21-801: 250ml	70mm	143mm	32s
Duran 21-801: 500ml	86mm	181mm	20s
Duran 21-801: 750ml	95mm	203mm	18s
Duran 21-801: 1000ml	101mm	230mm	10s
Duran 21-801: 2000ml	136mm	265mm	0

Vessel capacity without tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	50s
Duran 21-801: 150ml	62mm	110mm	45s
Duran 21-801: 250ml	70mm	143mm	32s
Duran 21-801: 500ml	86mm	181mm	20s
Duran 21-801: 750ml	95mm	203mm	18s
Duran 21-801: 1000ml	101mm	230mm	15s
Duran 21-801: 2000ml	136mm	265mm	8

s - staggered to fit maximum number