


Overall dimensions




Packed weight: 11kg
Product weight: 9kg

Packed dimensions;
A 660mm
B 490mm
C 330mm

Product dimensions;
A 545mm
B 380mm
C 225mm

Working Area



D	298 (5L) / 323 (12L) mm (at opening)
D	281 (5L) / 306 (12L) mm (at tray height)
E	148 (5L) / 298(12L) mm (at opening)
E	131 (5L) / 281(12L) mm (at tray height)
F	150mm (5L&12L)
G	132mm (5L&12L)
H	150mm (5L&12L)
I	191 mm (5L) & 241mm (12L)

Vessel capacity with tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	8 (5L) & 23s (12L)
Duran 21-801: 150ml	62mm	110mm	8 (5L) & 16 (12L)
Duran 21-801: 250ml	70mm	143mm	6s (5L) & 14s (12L)
Duran 21-801: 500ml	86mm	181mm	0 (5L) & 8s (12L)
Duran 21-801: 750ml	95mm	203mm	0 (5L) & 3s (12L)
Duran 21-801: 1000ml	101mm	230mm	0 (5L) & 2 (12L)
Duran 21-801: 2000ml	136mm	265mm	0 (5L) & 0 (12L)

Vessel capacity without tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	8 (5L) & 23s (12L)
Duran 21-801: 150ml	62mm	110mm	8 (5L) & 16 (12L)
Duran 21-801: 250ml	70mm	143mm	6s (5L) & 14s (12L)
Duran 21-801: 500ml	86mm	181mm	3 (5L) & 9 (12L)
Duran 21-801: 750ml	95mm	203mm	2 (5L) & 8s (12L)
Duran 21-801: 1000ml	101mm	230mm	2 (5L) & 5s (12L)
Duran 21-801: 2000ml	136mm	265mm	0 (5L) & 4 (12L)

s - staggered to fit maximum number