


Overall dimensions



Please note images are for representation purposes only


Packed weight: Tank: 3kg
Thermostat: 4kg
Total: 7kg

Product weight: Tank: 2.5kg
Thermostat: 2.8kg

Packed dimensions (h x w x d):
Thermostat: 310mm x 440mm x 200mm
Tank: 390mm x 250mm x 410mm

Product dimensions: A 235mm
B 325mm
C 180mm
D 385mm
E 172mm
F 120mm

Working area



G	155mm
J	120mm
H	150mm

Please note images are for representation purposes only

Vessel capacity with tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	4
Duran 21-801: 150ml	62mm	110mm	3s
Duran 21-801: 250ml	70mm	143mm	2
Duran 21-801: 500ml	86mm	181mm	0
Duran 21-801: 750ml	95mm	203mm	0
Duran 21-801: 1000ml	101mm	230mm	0
Duran 21-801: 2000ml	136mm	265mm	0

Vessel capacity without tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	4
Duran 21-801: 150ml	62mm	110mm	3s
Duran 21-801: 250ml	70mm	143mm	2
Duran 21-801: 500ml	86mm	181mm	1
Duran 21-801: 750ml	95mm	203mm	1
Duran 21-801: 1000ml	101mm	230mm	1
Duran 21-801: 2000ml	136mm	265mm	0

s - staggered to fit maximum number