


Overall dimensions



Please note images are for representation purposes only


Packed weight:	Tank: 4.5kg
	Thermostat: 4kg
	Total: 8.5kg
Product weight:	Tank: 3kg
	Thermostat: 2.8kg
	Total: 13.8kg

Packed dimensions (h x w x d):

Thermostat: 310mm x 440mm x 200mm
Tank: 330mm x 410mm x 410mm

Product dimensions:	A 185mm
	B 335mm
	C 215mm
	D 415mm
	E 172mm
	F 120mm

Working area



Please note images are for representation purposes only

G	150mm
H	130mm
J	150mm

Vessel capacity with tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	4
Duran 21-801: 150ml	62mm	110mm	3s
Duran 21-801: 250ml	70mm	143mm	2
Duran 21-801: 500ml	86mm	181mm	0
Duran 21-801: 750ml	95mm	203mm	0
Duran 21-801: 1000ml	101mm	230mm	0
Duran 21-801: 2000ml	136mm	265mm	0

Vessel capacity without tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	4
Duran 21-801: 150ml	62mm	110mm	3s
Duran 21-801: 250ml	70mm	143mm	2
Duran 21-801: 500ml	86mm	181mm	1
Duran 21-801: 750ml	95mm	203mm	1
Duran 21-801: 1000ml	101mm	230mm	1
Duran 21-801: 2000ml	136mm	265mm	0

s - staggered to fit maximum number