


## Overall dimensions



Please note images are for representation purposes only


<b>Packed weight:</b>	Tank: 4.5kg Thermostat: 4kg Total: 8.5kg
<b>Product weight:</b>	Tank: 3kg Thermostat: 2.8kg Total: 5.8kg

### Packed dimensions (h x w x d):

Thermostat: 310mm x 440mm x 200mm  
Tank: 330mm x 410mm x 410mm

<b>Product dimensions:</b>	A 185mm B 335mm C 215mm D 415mm E 172mm F 120mm
----------------------------	--

## Working area



<b>G</b>	<b>150mm</b>
<b>H</b>	<b>130mm</b>
<b>J</b>	<b>150mm</b>

Please note images are for representation purposes only

## Vessel capacity with tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	4
Duran 21-801: 150ml	62mm	110mm	3s
Duran 21-801: 250ml	70mm	143mm	2
Duran 21-801: 500ml	86mm	181mm	0
Duran 21-801: 750ml	95mm	203mm	0
Duran 21-801: 1000ml	101mm	230mm	0
Duran 21-801: 2000ml	136mm	265mm	0

## Vessel capacity without tank lid

	Diameter	Height	Max. no in bath
Duran 21-801: 100ml	56mm	105mm	4
Duran 21-801: 150ml	62mm	110mm	3s
Duran 21-801: 250ml	70mm	143mm	2
Duran 21-801: 500ml	86mm	181mm	1
Duran 21-801: 750ml	95mm	203mm	1
Duran 21-801: 1000ml	101mm	230mm	1
Duran 21-801: 2000ml	136mm	265mm	0

s - staggered to fit maximum number